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
Instructions for use

Please loosen the screws at the top and the bottom with the enclosed o-rings of the stirrer seal in order to introduce the glass or metal stirrer shaft. The shaft has to be dry, clean and free of oil or grease.

After inserting the shaft tighten both screws by hand and combine the mechanical seal with the glass part. A Viton-/Perfluor*-O-Ring is sealing automatically against the glass, when the screw thread is tightened.

Now you can connect the stirrer shaft with the stirrer motor by a stirrer coupling.

In order to prevent an accident, covering of rotating parts is recommended.

An elevated temperature of the housing of the mechanical seal is normal.

The mechanical seal is in principle appropriate for unlubricated operation.

Special operating conditions can make the admission of the sliding rings with a pressure-free foreign medium necessary (**Quench**). This can be brought in by the lateral drilling at the lower end of the mechanical seal housing. As a function of the respective requirement profile the following Quenchmediums can be used:

- **Liquid-Quench** (e.g. for cooling/lubrication of the sliding rings or for exclusion of air): If the temperature of the Quenchmediums is $< 60^{\circ}\text{C}$ e.g. water of low hardness (max. 5°dH) or VE Wasser/Glycerin, and/or Ethylenglycol (without inhibitors) 50/50 can be used.
At higher temperatures preferably light oil is to be used after arrangement.
- **Gas-Quench** (e.g. exclusion of air): e.g. nitrogen, dry air for anti-icing, etc..

The selection of the Quenchmedium depends also on other aspects like: Product compatibility, environmental compatibility, boiling point as high as possible, high specific thermal capacity, low viscosity.

A general maximum permissible RPM for non lubricated use could not be given, because it depends on various effects, like stirrer dimension, ambient temperature, experiment length, amount of wear, medium, pressure/vacuum**, etc.. This means the maximum RPM for dry operation depends on the experimental arrangement and must be tested.


If there are problems, after a longer period, like extreme mechanical seal temperatures or unregular noises please send us the whole mechanical seal for repair. Do not try to repair the mechanical seal yourself, because certain experience is required to void leakage or becoming tight.

Technical data

Parts in contact:	PEEK, Hastelloy, Viton or Perfluor*
Secondary shaft seal:	Viton or Perfluor*
Body:	PEEK, Aluminium
max. RPM:	$\sim 600 \text{ min}^{-1}$; short time (1-2h) $\sim 800 \text{ min}^{-1}$
temperature range.:	up to 200°C
shafts:	glass- or metalshaft, also coated, with dia. 8/10/14/16 mm

* various executions available

** The mechanical seal is for use under normal pressure or vacuum. By use under pressure a sealing up to 0,3-0,5 bar pressure can be reached.

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