


[Luer adapter for GL 14 - SON051121]

GL threads are used on glass apparatus in the laboratory. Adapters are required for the connection with smaller tubing diameters.

This adapter makes it possible to safely connect Luer connections with GL threaded fittings.

Due to the large number of different parts that can be connected via micro tubing, a high degree of variation is given and thus a large number of applications are possible. Since the "Luer" system is standardized, many parts can be combined with each other in any way.

Technical or chemical applications can be realized with this adapter. Connecting micro tubing to glass apparatus is made easy.



Item No.: SON051121

Consisting of:

- SON051121 - Luer adapter for GL 14

Necessary accessories:

- 132922705 - Cap GL 14

Alternative accessories:

- 106720109 - Tubing 1,5 x2,7 mm, length 1,5 m
- 19421100*(4-6) - Cannulas (selection, depending on application)

Alternative Adapters:

- SON051131 - Luer adapter for GL 18