

Pipette Tips, Eppendorf

The pipette tips for the Varipette® 4720 are tailored to different vessels. The Varitip P is designed for aspirating 1 ml to 10 ml, for example, from beaker. It pipettes according to the positive-displacement principle. The Varitip S for 2.5 ml to 10 ml forms a system with the Maxitip and is used for aspirating liquid from tall, narrow-neck vessels.

- Varitip P for aspirating 1 ml to 10 ml from beakers – positive-displacement principle
- Varitip S for aspirating 2.5 ml to 10 ml from tall, narrow-neck vessels. System consists of a dispensing part and the Maxitip - air-displacement principle
- No carryover as Maxitip exchanged
- Valve for Maxitip ensures drip-free dispensing of liquids with a high vapor pressure
- Graduated Maxitip G for compensation of nominal value of Varipette

Varitips P

Article	Order No.
Starter kit	14 69 47208
100 Maxitips, 10 dispensing parts, 10 valves, 2 Maxitips G	
Varitips P	14 69 47202
For aspirating from smaller vessels, set of 100	

Varitips S

For aspirating from narrow-neck vessels and volumetric flasks, consisting of:

Article	Order No.	PU pc/Pack
Dosing part for Varitips S	14 69 47205	30
Maxitip for Varitips S	14 69 47206	200
Maxitip-G, graduated tip for Varitips S	14 69 47209	5
Valve for Varitips S	14 69 47207	100

❗ Available only in original packaging units. No fractions possible.

