

Bottle-Top Burette, Eppendorf

The innovative Top Buret from Eppendorf sets new standards for manual titration. A new, pulse, free dispensing technique allows the desired volumes to be dispensed continuously, thus eliminating the need for new liquid aspiration each time. This makes titration simpler and safer, with precision values well within required limits.

The medium is dispensed by simply rotating the wheels, with a dispensing rate of 2,500 μl per rotation for the model M and 5,000 μl per rotation for the model H.

- Continuous, pulse-free dispensing technique for rapid titration
- Patented safety valve for deairing without loss of reagent
- Voltage supply with long-life primary cells and charging display
- Dispensing range from 0.01 ml to 999.9 ml
- Removable electronics unit
- Modular, service-friendly construction
- Simple calibration program
- Discharge tube rotatable by 360°
- Adjustable discharge tube horizontal: 142 to 220 mm, vertical: 10 to 200 mm
- Telescopic aspiration tube for automatically setting the length from 210 mm to 370 mm
- Connection for dry tube

Top Buret with safety valve, telescopic aspiration tube, adjustable discharge tube, three adapters for 40/38/32 mm bottle threads, 2 $\times$ 1.5 V microbatteries.

Article	Order No.
Eppendorf Top Buret M, 25 ml	14 69 50017
Eppendorf Top Buret H, 50 ml	14 69 50025
Dry tube	14 69 51000

