

Separatory Funnels, NALGENE®, IDL

Separatory funnels made of polypropylene with stopcock and plug made of PTFE TFE and polypropylene screw closure. The phase interface is visible down to the stopcock. Leakproof screw closure. Resists virtually all chemicals, including hydrofluoric acid. Stopcock and plug are easily removed for cleaning, at which point the funnels can be sterilized in an autoclave.

Nominal volume ml	Cap Ø mm	Stem length mm	Order No.
125	28	65	23 87 43001
250	33	65	23 87 43002
500	43	65	23 87 43003
1000	53	65	23 87 43004

N.B.: Do not autoclave when closed! (Make sure the stopcock or closure is open).

