

SnapIt Ampoule Opener, neoLab®

SnapIt – A simple safer way to open glass ampoules. Protect yourself from the risk of contamination and sharp injury. Insert the ampoule lid into the SnapIt, snap open the ampoule and eject the ampoule lid directly into a sharps bin. The SnapIt Ampoule opener is easy to clean and can be autoclaved if required. The flat edge of the Trolley edition prevents an inadvertent rolling.

Article	Colour	for volume ml	Order No.
SnapIt Personal, regular	blue	1 – 15	10 01 20040
SnapIt Trolley, regular	black	1 – 15	10 01 20041
SnapIt Trolley, large	red	5 – 25	10 01 20042

