

Agate Balls for Laboratory Mills, IDL



Ball mixes for agate milling vessels:

ml	Ball dia.	3 mm	5 mm	7 mm	10 mm	15 mm	20 mm	25 mm	30 mm	40 mm
50	qty	400	80	35	10	4	2	–	–	–
602	qty	2000	400	150	40	18	8	4	3	–
602	qty	–	800	400	100	40	18	10	7	4

Ø approx. mm	Order No.	Ø approx. mm	Order No.
3	21 00 60403	20	21 00 60420
5	21 00 60405	25	21 00 60425
7	21 00 60407	30	21 00 60430
10	21 00 60410	40	21 00 60440
15	21 00 60415		

Other sizes on request.

