

Restrainer for Magnetic Stirring Bars, neoLab®

Easily retrieve magnetic stirring bars from vessels!

This restrainer **stops the loss of magnetic stirring bars** headed down lab drains and eliminates the splashing hazard when pouring out of mixing containers. This easy to use tool fits comfortably in one hand and when placed next to the outside of the container, readily attracts and holds stirring bars in place. It contains a powerful magnet. As containers are emptied, the stirring bar is retained and can be safely removed from the vessel.

- easy to use
- for all magnetic stirring bars up to 80 mm
- easy to clean

Material: polystyrene

Article

Restrainer for magnetic stirring bars

Order No.

10 01 61179

