

Colorization Bank, Landgraf

Stainless steel Colorization Bank with tipping device for staining of slides for laboratory work. The Colorization Bank has 2 independently operated tipping devices for slides. The tipping device can be set for staining and for washing the slides. The waste water is collected in a sump and can be drained away into a sink.

The front section of the bank can be folded back, so that the slides can be heated up with a Bunsen burner (e.g. for TB staining).

Comprising: Colorization Bank with collection tank and stand

Article	Working W mm	Order No.
Colorization bank	330	88 67 30330
Colorization bank	450	88 67 30450
Colorization bank	570	88 67 30570

