

Counters, Hecht-ASSISTENT CE

The AC-12 and AC-15 electronic memory counters with 8, 12 or 15 keys simplify blood count analysis. Counting is stopped automatically upon reaching a full hundred. Counting beyond a hundred is possible to a maximum of 1000.

Counting can be stopped at any number required. The stored results can be displayed as absolute or percentage values, with switching between the two display modes as often as required.

Briefly pressing the Clear key has no effect, thus ensuring that no data is deleted unintentionally. The memory is not cleared until the key is kept depressed for approx. 1 to 2 seconds.

The last entry in a sub-memory can be deleted from the sub-memory and the total memory using the "error corr".

Technical data:

Dimensions: 19,5 x 15 x 3 cm
Power supply: 220 V, 50 Hz, 3 W (COUNTER AC-8 and -12)
220 V, 50 Hz, 6 W (COUNTER AC-15)

Article	Order No.
① Counter AC-8	87 18 10008
② Counter AC-12	87 18 10012
Counter AC-12 with built-in rechargeable battery pack	87 18 10013
Counter AC-15	87 18 10015
numeric keys with illuminated display, 4 defined programs and a free counting program	
Counter AC-15 with PC- interface	87 18 10016
Microsoft Windows 2000 or XP Home/Professional	

