

UV PCR Workstation, LAT

Shortwave UV irradiation and antimicrobial metal surfaces for a working environment securely protected against contamination.

When the lid is opened, the UV light is switched off automatically. Two power outlets in the interior are provided for the connection of small devices, two shelves for storage of small objects. With the integrated timer function, the device can be sterilized between experiments. White light can be enabled and disabled separately. The UV2 model has additional UV light sources in the air recirculation.

Article	Dimensions (W×H×D) mm	Order No.
UV PCR Cabinet	544×610×729	10 15 27121
UV2 PCR Cabinet	544×610×729	10 15 27122

Also available as a larger PCR Workstation with four power outlets!

