

Centrifuge Bottles, Thermo Scientific Nalgene, IDL



WWRJN117

thermo scientific

Autoclavable centrifuge bottles made of polypropylene copolymer with polypropylene screw closure. Translucent bottles have excellent chemical resistance. Strong and rigid. Sealing cap available for 1000-ml capacity bottles. When using the 250-ml bottle in the Sorvall GSA rotor, use the NALGENE centrifuge bottle adapter. Order No. 23 87 31202 is designed for use with IEC rotors. Before autoclaving, set closure on top of the bottle without engaging the threads. The 2000-ml bottle is also available as a sterile version; for both variants, the support adapter (Art. No. 75007686) is required, see picture. For proper performance these bottles must be filled to at least 80% of the total capacity, except for Order No. 23 87 31203, which must NOT be filled to more than 75% of total capacity.

Volume ml	load capacity up to g	Order No.	PU pc/Pack
250	13200	23 87 31201	4
500	4800 (for IEC rotors)	23 87 31202	4
500	4800	23 87 31203	4
1000	7100	23 87 31204	4
1000	7100 (for IEC rotors)	23 87 31205	4
2000	7333 (for H-12000 rotors)	10 87 31206	6
2000 (steril)	7333 (for H-12000 rotors)	10 87 31207	24

❗ Available only in original packaging units. No fractions possible.

N.B.: An adapter is essential for Sorvall GSA rotors (see Order No. 23 87 31250).

