

Heat-Transfer Fluids LABOTHERMOL®, neoLab®

For oil baths and melting point measuring equipment, for gentle heating, for prevention of overheating, for heat transfer in the lab and workshop, according to requirements.

Article

Order No.

- ① **LABOTHERMOL® O, heating bath liquid on a mineral oil basis** 35 01 16060

dark colored, high resistance to oxidation, canisters of 5 l

Operating range: -10 to +280°C

Viscosity at 40°C: 43.5 mm²s⁻¹ (20 cSt.)

Article

Order No.

- ② **LABOTHERMOL® S, silicone oil heating bath liquid** 35 01 16070

colorless and hence also ideal for melting point measuring equipment. With stabilizing additive, bottles of 1000 ml

Operating range: 0 to +150°C

Viscosity at 25°C: 1000 cSt.

Article

Order No.

- ③ **LABOTHERMOL® K, safety cold and heat transfer fluid on an alcohol basis** 35 01 16065

nearly colorless, water-soluble and hence washable, canisters of 5 l

Operating range: -20 to +150°C

Flashpoint: >+255°C acc. to DIN 51758



①



③



②

