

Eye Wash Station, WINLAB®

For optimum first-aid in case of eye accidents.

For effective eye washing in just 1 to 2 steps!

Complete with a 200 ml sterile pH-NEUTRAL bottle as an effective phosphate buffer in case of accidents involving acids, and 500 ml 0.9% sodium chloride solution for further washing.



WW5PT060

WINLAB®

- | Article | Order No. |
|--|--------------------|
| ① Emergency eye wash station on a wall holder with mirror
Thanks to the wall holder with mirror, both bottles are INSTANTLY TO HAND! | 36 43 60008 |
| ② DUO emergency eye wash station with wall holder
Eye wash station DUO for simultaneous washing of both eyes. With a pH-neutral DUO eye wash bottle containing 500 ml of sterile phosphate buffer solution (4.9%) and a DUO eye wash bottle with 1,000 ml sodium chloride solution (0.9%). With wall mount, pictogram, and mirror. Eye emergency care plan is comprised in the delivery. Suitable in particular for use at workstations where the eyes may come into contact with acids, bases, and foreign particles. | 10 03 60014 |
| ③ Emergency eye wash station in a dust-proof wall box for flipping open
With mirror and separate adhesive sign with directions for use. Wall box 265×227×110 mm (H×W×D), sign 290×90 mm. | 36 43 60009 |
| ④ DUO emergency eye wash station in a dust-proof wall box
Eye wash station DUO for simultaneous washing of both eyes, supplied in a dust-proof wall box. With a pH-neutral DUO eye wash bottle containing 500 ml of sterile phosphate buffer solution (4.9%) and a DUO eye wash bottle with 1,000 ml sodium chloride solution (0.9%). With wall mount, pictogram, and mirror. Suitable in particular for use at workstations where the eyes may come into contact with acids, bases, and foreign particles. | 10 03 60016 |

