

Cleaning Detergent LABOSOL® S, neoLab®

Powerful concentrated cleaner for dirty laboratory glass devices, porcelain, instruments and materials made of plastic and rubber.

Thorough, self-acting cleaner. Does not attack the material being cleaned or the skin. Also suitable for use in ultrasonic cleaners.

Use a 0.5 – 2% solution with warm water depending on how dirty the items are. LABOSOL® S contains no chlorine or caustic potash, is bio-degradable and can be used as a substitute for a chromic-sulfuric acid mixture.

Article	Volume l	Order No.
① Labosol® S, 5 liter	5	35 01 16090
Labosol® S, 1 liter	1	35 01 16091

