

Nitrogen Vacuum Workplaces, Rettberg



WWQVD181

rettberg.biz
GLASAPPARATEBAU • LABORATORY GLASSWARE

Complete unit (installation width approx. 115 cm), consisting of:

- Shutoff valve, 0 to 15 mm, removable, with KF16 for suction-side connection to double cooling trap with additional 0 to 8 mm valve and 10 mm olive for vacuum or ventilation
- Vacuum double cooling trap directly connected, KF16 connections at the side, bottom elements KF40 with approx. 350 ml volume. Length of bottom elements approx. 210 mm
- N2/vacuum manifold KF16 with 5 two-way valves 0 to 8 mm, connection for N2 or vacuum via KF16 or 10 mm olive



Article	Order No.
① Complete workstation	10 70 00222
N2/vacuum interchangeable rack	10 70 00223
Vacuum double cooling trap	10 70 00224
Vacuum shut-off valve	10 70 00225

