

Dry Ice Machine SnowPack®, BÜRKLE



WWR4Q053



SnowPack® enables you to obtain a -79 °C dry ice tablet weighing 50 or 250 g within a few moments, without complicated technology and without electricity!

SnowPack® is used in hospitals, universities, schools, institutes, laboratories, in the chemical industry and in materials testing facilities.

Examples of use: lowering the temperature in cooling baths or cold traps, pre-cooling tissue samples, ensuring safe transport or embedding biological substances, freezing blood and tissue samples, studies in chemistry, physics and materials technology, etc.. The SnowPack® dry ice machine can be installed on any commercial CO₂-bottle with a standpipe (without pressure regulator) (thread W 21.80 × 1/14" acc. to DIN 477 Nr. 6).

Article	Tablets weight approx. g	Tablet dia. × H mm	Order No.
SnowPack® 50	50	50 × 30	10 04 59853
SnowPack® 250	250	80 × 60	10 04 59852

