

Filter Papers for Sugar Analysis, Hahnemühle

Recommended for the clarifying filtration of beet extracts and juices prior to analysis.

- High filtration speed combined with high retention of particles
- Two surfaces are available: smooth or creped

Technical data:

Grade	Properties	Filtration Herzberg s	Thickness mm
3459	fast, creped	110	0.3
3002	medium fast, smooth	150	0.14

Grade	Ø mm	Weight g/m ²	Order No.	PU pc/Pack
3002	200	69	10 15 93002	1000
3459	230	74	10 15 93459	100

Other sizes on request.



WW6PB017



Hahnemühle

