

Filter Head with PP-Funnel and 2 FKM Seals, DURAN®, DURAN Group / IDL



DURAN® filter heads with PP funnel and 2 FKM seals

Filter disks must always be ordered separately.

Important: The filter disk has to be placed between 2 seals.

Easy, safe removal of the filtered material. Longer useful life of the filter disk because the filtered material is not scraped off. Accessible from both sides for easy cleaning.

Volume ml	for plate Ø mm	Thread mm	Order No.
30	24	28	11 24 72024
250	50	54	11 24 72050
1000	90	95	11 24 72090

① Single parts: DURAN® filter head with thread

Volume ml	Thread mm	Order No.
30	28	11 24 72124
250	54	11 24 72150
1000	95	11 24 72190

② FKM seals (fluorocautouc)

for plate Ø mm	Order No.	PU pc/Pack
24	11 29 22024	10
50	11 29 22050	10
90	11 29 22090	10

❗ Available only in original packaging units. No fractions possible.

③ Filter disks with glass rim

Plate Ø mm	Porosity	Order No.	Plate Ø mm	Porosity	Order No.
24	1	11 25 20241	50	3	11 25 20503
24	2	11 25 20242	50	4	11 25 20504
24	3	11 25 20243	90	1	11 25 20901
24	4	11 25 20244	90	2	11 25 20902
50	1	11 25 20501	90	3	11 25 20903
50	2	11 25 20502	90	4	11 25 20904

④ PP funnel, autoclavable up to +140 °C

Thread mm	Stem O.D. mm	Order No.
28	10	11 29 22124
54	12	11 29 22150
95	18	11 29 22190

⑤ DURAN® slotted filter disk (holder for membrane and paper filters)

Fits for order No. 112472050

for plate Ø mm	Order No.	PU pc/Pack
48	11 21 34031	10

❗ Available only in original packaging units. No fractions possible.

