

## Inoculating Loops, Stands and Holders, neoLab®

### Inoculating holders

Kolle needle holders for inoculation loops, color-coded.

| Article                        | Colour | Overall length approx.<br>mm | Order No.          |
|--------------------------------|--------|------------------------------|--------------------|
| Inoculating loop holder white  |        | 245                          | <b>33 01 12120</b> |
| Inoculating loop holder blue   |        | 245                          | <b>33 01 12121</b> |
| Inoculating loop holder black  |        | 245                          | <b>33 01 12122</b> |
| Inoculating loop holder red    |        | 245                          | <b>33 01 12123</b> |
| Inoculating loop holder yellow |        | 245                          | <b>33 01 12124</b> |

### ① Inoculating stands, rectangular

| Article                | Material  | for qty inoc. loops | L x W x H<br>mm | Order No.          |
|------------------------|-----------|---------------------|-----------------|--------------------|
| Inoculating loop stand | Aluminium | 4                   | 140 x 40 x 60   | <b>33 01 12118</b> |
| Inoculating loop stand | Wood      | 4                   | 140 x 40 x 60   | <b>33 01 12125</b> |

### ② Inoculating stands, round

Including space for inoculating loop storage tube.

| Article           | Material  | for qty inoc. loops | Ø<br>mm | Height<br>mm | Order No.          |
|-------------------|-----------|---------------------|---------|--------------|--------------------|
| Inoculating stand | Aluminium | 6                   | 100     | 60           | <b>33 01 12119</b> |
| Inoculating stand | PP        | 6                   | 100     | 60           | <b>33 01 12129</b> |

### Inoculating loops made of special stainless steel

Suitable for repeated annealing.

| Article          | Material        | Ø loop<br>mm | L<br>mm | Order No.          |
|------------------|-----------------|--------------|---------|--------------------|
| Inoculating loop | Stainless steel | 1            | 50      | <b>33 01 12115</b> |
| Inoculating loop | Stainless steel | 3            | 50      | <b>33 01 12116</b> |
| Inoculating loop | Stainless steel | 5            | 50      | <b>33 01 12117</b> |

### Inoculating loops made of platinum/iridium wire

| Article          | Material         | Ø wire<br>mm | Ø loop<br>mm | L<br>mm | Order No.          |
|------------------|------------------|--------------|--------------|---------|--------------------|
| Inoculating loop | Platinum iridium | 0.6          | 1            | 50      | <b>33 01 12111</b> |
| Inoculating loop | Platinum iridium | 0.6          | 3            | 50      | <b>33 01 12112</b> |
| Inoculating loop | Platinum iridium | 0.6          | 5            | 50      | <b>33 01 12113</b> |
| Inoculating loop | Platinum iridium | 1.0          | 3            | 50      | <b>33 01 12127</b> |
| Inoculating loop | Platinum iridium | 1.0          | 5            | 50      | <b>33 01 12128</b> |

