

## Distilling Systems, behrotest®, Behr

### behrotest® distilling apparatuses for determination of the phenol index by the DIN/DEV 38409-H16-2 method

DIN/DEV requires the distillation of a large amount of liquid within a short time for the determination of the phenol index.

The behrotest® WE 5 five-output steam generator can supply a maximum of 5 samples simultaneously for fast and gentle distillation. Steam input into the samples guarantees uniform distillation and high reproducibility of retrieval rates over the entire range.

The five-output linear steam generator with central water supply and level control was specially developed for this application.

| Article                                                                          | Order No.          |
|----------------------------------------------------------------------------------|--------------------|
| <b>Complete distilling apparatus for the determination of phenols</b><br>5 place | <b>10 48 91229</b> |

Supplied without the illustrated accessories

