

Delbrück Capillaries end-to-end DURAN, Hirschmann

End-to-end capillaries according to the Prof. Delbrück method, Na-heparinized, with declaration of conformity acc. to DIN 12600, with batch number on packaging, with CE mark, acc. to IVD regulations and ISO 7550. Unit pack 10 × 100 pieces.

Volume µl	Variation coefficient %	Accuracy %	Length mm	Order No.	PU pc/Pack
20	1.0	0.5	30	19 37 20301	1000
40	1.0	0.5	30	19 37 20302	1000
50	1.0	0.5	30	19 37 20303	1000



HIRSCHMANN®

webcode  WW4ZV051

