

Cooling glass vessels

For use with HD 2000 / HD 3000 / HD 4000

Order no. KG 3 - 536 / KG 5 - 481



The cooling glass vessels KG 3 and KG 5 with cooling jacket are especially suited for processing temperature-sensitive samples.

The coolant passing the cooling jacket enables temperature constancy during sonication.

	KG 3	KG 5
probe Ø:	13 mm	25 mm
HD 2070 / HD 3100 / HD 4050 / HD 4100	to 15 ml	-
HD 2200 / HD 3200 / HD 4200	to 15 ml	to 70 ml

Supply

1 cooling vessel made of glass (KG 3 or KG 5)

	KG 3	KG 5
Dimensions, interior	Ø 20 x 60 mm	Ø 35 x 100 mm
Diameter, exterior	Ø 40 mm	Ø 55 mm
Olive fittings	Ø 10 mm	Ø 12 mm
Height, total	75 mm	115 mm

The measure data are only approximate values.

When positioning and adjusting the immersing depth, avoid contact of probes with the vessel - danger of breakage!

