

## Benchtop Centrifuge ROTOFIX 32 A, HETTICH

### Benchtop centrifuge ROTOFIX 32 A

Microprocessor-controlled bench-top centrifuge with a large and diverse range of accessories. With frequency drive, lid catch and lock, impulse key for fast spin, one-hand closing, automatic rotor recognition and imbalance cut-out. Easy input of speed and time on the membrane keypad. Display of ACTUAL status during operation.

#### Technical data:

|                              |                         |
|------------------------------|-------------------------|
| External dimensions (WxDxH): | 366 × 430 × 257 mm      |
| Weight:                      | approx. 23 kg           |
| Power supply:                | 208 – 240 V, 50/60 Hz   |
| Connected load:              | 300 VA                  |
| Max. speed:                  | 6000 min <sup>-1</sup>  |
| Max. RCF:                    | 4226                    |
| Max. capacity:               | 4 × 100 ml / 32 × 15 ml |
| Emission/Immunity:           | EN/IEC 61326-1, class B |



①

| Article                                           | Order No.   |
|---------------------------------------------------|-------------|
| ① ROTOFIX 32 A, benchtop centrifuge without rotor | 10 56 11206 |

#### Accessories

| Article                                                                                                   | Order No.   |
|-----------------------------------------------------------------------------------------------------------|-------------|
| ② <b>Swig-out rotor, without suspension gear</b><br>for max. 4 × 100 ml or 4 slide samples, max. RCF 2665 | 55 56 11624 |
| <b>Round bucket, alloy, 50 ml, fits in rotor 1324</b>                                                     | 10 56 11398 |
| <b>Round bucket, alloy, 100 ml, fits in rotor 1324</b>                                                    | 10 56 11390 |
| <b>Round bucket, alloy, 100 ml, fits in rotor 1624</b>                                                    | 55 56 11381 |
| <b>Lid for bucket 1381, 1430 and 1390</b>                                                                 | 55 56 11382 |
| ③ <b>Slotted suspension gear, fits in rotor 1624</b>                                                      | 55 56 11742 |
| <b>Swing-bucket rotor for max. 8 × 15 ml, without suspension gear</b>                                     | 55 56 11611 |
| ④ <b>Swing-bucket rotor for max. 8 × 50 ml conical tubes</b><br>with screw cap, with suspension gear      | 55 56 11617 |
| ⑤ <b>Swing-bucket rotor for max. 6 × 50 ml conical tubes</b><br>with screw cap, with suspension gear      | 55 56 11619 |
| <b>Swing-bucket rotor for max. 12 × 15 ml, without suspension gear</b>                                    | 55 56 11628 |
| <b>Swing-bucket rotor for max. 4 × 100 ml, without suspension gear</b>                                    | 10 56 11324 |
| ⑥ <b>Fixed-angle rotor for max. 6 × 85 ml</b>                                                             | 10 56 11620 |
| <b>Fixed-angle rotor for max. 12 × 15 ml, max. RCF 4146</b>                                               | 55 56 11613 |
| <b>Fixed-angle rotor for 32 × 15 ml</b>                                                                   | 10 56 11418 |
| ⑦ <b>Zyto rotor for 6 microscope slide samples, without suspension gear</b>                               | 55 56 11626 |
| <b>Zyto rotor for 8 microscope slide samples, with suspension gear</b>                                    | 55 56 11648 |

More accessories on request.



②



③



④



⑤



⑥



⑦