

Sample Splitters, HAVER & BOECKER

Representative test samples can be easily and reliably produced with HAVER sample splitters

① Sample splitters HAVER RT 6,3 / RT 12,5

1 rack: 300 × 250 × 270 mm, made of lacquered sheet steel
1 divider: openings depending on models, stainless steel
3 receivers: contents 1.75 litres, made of tinned sheet

Article	Version	Order No.
RT 6.3	12 divisions 6.3 mm / ¼ inch	10 16 90292
RT 12.5	18 divisions 12.5 mm / ½ inch	10 16 90291

② Sample splitters HAVER RT 25 / RT 37.5 / RT 50 / RT 75

1 rack: 620 × 260 × 420 mm, galvanized sheet steel
1 divider: openings depending on models, galvanized sheet steel
3 receivers: contents 8 litres, galvanized sheet steel

Article	Version	Order No.
RT 25	16 divisions 25 mm / 1 inch	10 16 90290
RT 37.5	18 divisions 37.5 mm / 1½ inch	10 16 90288
RT 50	8 divisions 50 mm / 2 inch	10 16 90287
RT 75	6 divisions 75 mm / 3 inch	10 16 90286

③ Rotating Sample Reducer HAVER RPT

Splitting channels: depending on models 8, 10 or 30
Splitting steps: 3000 per minute
Max. feed volume: 4000 ml (RPT 1:8), 2500 ml (RPT 1:10) or 300 ml (RPT 1:30) per run
Max. particle size: 10 mm (RPT 1:8 and RPT 1:10) or 2.50 mm (RPT 1:30)

Article	Order No.
HAVER RPT 1:8 (8 glass sample bottles)	10 16 90991
HAVER RPT 1:10 (10 glass sample bottles)	10 16 90992
HAVER RPT 1:30 (3 glass sample bottles)	10 16 90994

