

Erlenmeyer Flasks, DURAN®, DURAN Group / IDL

Erlenmeyer flasks with DIN threads and graduation, without illustrated caps

With Retrace Code (Batch Certificate of Performance via Internet).

Volume ml	Thread GL	Flask max. outer Ø mm	Height mm	Order No.	PU pc/Pack
100	25	64	105	11 21 80324	10
250	32	85	145	11 21 80336	10
500	32	105	175	11 21 80344	10
1000	32	131	220	11 21 80354	10

Caps see page XXX.

Important: Autoclave these flasks only with a loosely screwed cap (do not close tightly).

Erlenmeyer flasks with DIN threads and graduation, with red PBT-caps

With Retrace Code (Batch Certificate of Performance via Internet).

Volume ml	Thread GL	Flask max. outer Ø mm	Height mm	Order No.	PU pc/Pack
100	25	64	109	10 21 15324	10
250	32	85	149	10 21 15336	10
500	32	105	180	10 21 15344	10
1000	32	131	225	10 21 15354	10

