

Cooling Traps, Rettberg

Cooling Traps

Lower section removable, complete unit, suitable for use in vacuum conditions, 3.3 borosilicate glass, capacity approx 500 ml (condensate liquid approx 200 ml).

Lower section with small flange, 50 mm, seal and clamp.

Upper section with 2 small lateral flanges 16 mm (thick-walled), for fitting to high vacuum leads, metal flange components and similar.

Article	Order No.
Cooling trap, 2-piece	13 40 51144
Cooling trap 2-part	13 40 51244

with an additional small flange DN 16 (standard flange) for fitting a Pirani probe

Supplied without the stand shown in the illustration.

Accessories

Article	Order No.
Set of connecting pieces (2 pieces)	13 40 51316
small flange KF 16 mm (thick wall), on 16 mm hose shaft with clamp and seal (for installing a hose line)	
Clamp-ring holder	13 40 51410
with rod for fastening the cooling trap to a stand	
Dummy cap for closing the Pirani connection	13 40 51416
with small flange DN 16 mm, clamp and seal	

