

Cooling Traps, Rettberg

Cooling Traps

Lower section removable, complete unit, suitable for use in vacuum conditions, borosilicate glass 3.3, capacity approx 1000 ml (liquid condensate volume approx 500 ml)

Lower section with flange DN 60, seal and clamp

Upper section with 2 small lateral flanges DN 25 mm (thick-walled), for fitting to high vacuum leads, metal flange components and similar.

Article	Order No.
Cooling trap 2-piece	13 40 51154
Cooling trap 2-part	13 40 51254

with an additional small flange DN 16 (standard flange) for fitting a Pirani probe

Supplied without the stand shown in the illustration.

Accessories

Article	Order No.
Set of connecting pieces (2 pieces)	13 40 51425

Small flange KF 25 mm (thick wall), on 16 mm hose nozzle with clamp and seal (for installing a hose line)

Clamp-ring holder	13 40 51420
--------------------------	--------------------

with rod for fastening the cooling trap on a stand

Dummy cap for closing the Pirani connection	13 40 51416
--	--------------------

with small flange DN 16 mm, clamp and seal

