

Eye Wash Bottles, WINLAB®

Ready within seconds to protect your eyes!

This fast, flexible and simple eye wash is a sterile solution (0.9% salt solution) which prevents dirt and dust particles or micro-organisms from penetrating further into the eye. The eyewash bottles have an ergonomically shaped eye cover to keep the eyelids open during washing so that the eyewash gets directly to the eye. Unopened bottles keep for up to three years. No maintenance required.



WINLAB®

Article	Order No.
① Small bottle, 200 ml For personal protection kits.	36 43 60001
④ DUO Eye wash bottle, 1000 ml with 1,000 ml sodium chloride solution, for simultaneous washing of both eyes	10 03 60010
② Eye wash station, 500 ml With 1 bottle, a practical wall holder and pictures giving directions for use. Approx. 3 minutes rinsing time.	36 43 60003
⑥ Storage box for eye wash bottles WINLAB Midi Box (236×225×125 mm) for storing 200 ml and 500 ml eye wash bottles, also suitable for 500 ml DUO eye wash bottle	10 03 60011

Neutralizing eye rinse solution

For protection in case of accidents with acids and lyes. The eyewash bottle is filled with 4.9% phosphate buffer solution which neutralizes the harmful substances in acids and lyes and, if the eye is washed sufficiently, returns the pH level in the area of the eye to the normal pH of 7.3.

Article	Order No.
③ Small bottle, 200 ml For personal protection kits.	36 43 60005
⑤ DUO Eye wash bottle, 500 ml pH neutral, for simultaneous washing of both eyes	10 03 60012

