

Draining Racks, BÜRKLE-LaboPlast®

Plastic draining racks (PVC), largely chemically resistant, breakproof, used in the lab by the thousand and adapted to individual requirements.

The sturdy rods can be inserted in the holes as required and clicked in place. Using a combination of three different lengths of rod (60 mm, 100 mm and 150 mm) the rack can be adapted to take any item. The rods can be removed and repositioned at any time. With a removable draining channel plus a connection for a draining tube (to direct the water into the sink).

Available in 3 different sizes.

Supplied in strong cardboard boxes, complete with peel-off foil to protect from dust and scratching, rods and dowels and wall mounting.



WWR4Q022



①

① Draining racks with channel complete with rods

Length cm	Width cm	Qty rods 60 mm	Qty rods 100 mm	Qty rods 150 mm	Order No.
40	40	5	26	2	24 05 64040
50	50	10	39	4	24 05 64050
60	60	15	56	6	24 05 64060

Rods, replaceable/repositionable

Length mm	Rod Ø mm	Order No.
60	6	24 05 64460
100	10	24 05 64410
150 ¹⁾	12	24 05 64415

¹⁾ Rods only usable together with 100 mm.

② Draining rack, benchtop, with collecting pan complete with rods

The draining rack is easy and flexible to set down on a table or working surface. Two large supports with wide bases ensure sufficient stability. Since both sides of the draining board can be used, it provides twice the hanging space. The collecting pan can simply be removed for emptying.

Length cm	Width cm	Qty rods 60 mm	Qty rods 100 mm	Qty rods 150 mm	Order No.
40	40	5	26	2	10 04 53041
50	50	10	39	4	10 04 53051
60	60	15	56	6	10 04 53061



②

