

Cooling Traps, Rettberg

Cooling Traps

Lower section removable, complete unit, suitable for use in vacuum conditions, borosilicate glass 3.3, capacity approx 1000 ml (liquid condensate volume approx 500 ml)

Lower section with flange DN 60, seal and clamp

Upper section with 2 small lateral flanges DN 25 mm (thick-walled), for fitting to high vacuum leads, metal flange components and similar.

Article	Order No.
Cooling trap 2-piece	13 40 51154
Cooling trap 2-part	13 40 51254

with an additional small flange DN 16 (standard flange) for fitting a Pirani probe

Supplied without the stand shown in the illustration.

Accessories

Article	Order No.
Clamp-ring holder	13 40 51420

with rod for fastening the cooling trap on a stand



rettberg.biz
GLASAPPARATEBAU • LABORATORY GLASSWARE

Webcode  **WWQVD170**

