

Pipette Tips, Eppendorf



eppendorf

Features and Benefits of the entire tip generation:

- New tip design enables better access to a wider variety of tips and wells
- Ergonomic: New geometry requires less force for attaching and ejecting the tips
- Soft, flexible tip cone results in better pipette tip seal
- Available in three purity levels (Eppendorf Quality/PCR clean/biopur)

① Pipette tips epT.I.P.S.® Standard

- Original, high-quality Eppendorf pipette tips at an attractive price
- Packed in re-sealable bags
- Available in all tip sizes from 10 µl to 10 ml
- 200 µl, 300 µl and 1000 µl tips are color coded, yellow and blue, respectively

Volume	Length mm	Order No.	PU pc/Pack
0.1 – 10 µl	34	14 69 00811	1000
0.1 – 20 µl	40	14 69 00838	1000
0.5 – 20 µl	46	14 69 00854	1000
2 – 200 µl (2 – 100)	53	10 69 90889	1000
20 – 300 µl	55	10 69 30900	1000
50 – 1000 µl	71	10 69 30927	1000
50 – 1250 µl	76	14 69 00935	1000
50 – 1250 µl	103	10 69 00730	1000
250 – 2500 µl	115	14 69 00951	500
100 – 5000 µl	120	14 69 00978	500
0.2 – 5 ml	175	10 69 30650	300
1 – 10 ml	165	14 69 00765	200
1 – 10 ml	243	14 69 00781	200

❗ Available only in original packaging units. No fractions possible.

② Pipette tips epT.I.P.S.® Box

- epTIPS Box: Working box full of tips
- Guard lid on refill trays enables contamination-free transfer of trays to the working box (epT.I.P.S. Box)
- System optimized for use with multi-channel pipettes
- Far less waste compared to disposable boxes. All components 100% recyclable

Volume µl	Length mm	Order No.	PU pc/Box
0.1 – 10	34	10 69 73010	96
0.1 – 20	40	10 69 73037	96
0.5 – 20	46	10 69 73053	96
2 – 200 (2 – 100)	53	10 69 03070	96
20 – 300	55	10 69 73096	96
50 – 1000	71	10 69 73118	96
50 – 1250	76	10 69 73134	96
50 – 1250	103	10 69 73657	96
500 – 2500	115	10 69 73150	48
100 – 5000	120	10 69 73177	24

❗ Available only in original packaging units. No fractions possible.

③ Pipette tips epT.I.P.S.® Reloads and epT.I.P.S.® Reloads PCR clean

- epTIPS Reloads: Refill trays for the working box
- Guard lid on refill trays enables contamination-free transfer of trays to the working box (epT.I.P.S. Box)
- Refill trays (epT.I.P.S. Reloads) can be autoclaved in advance
- System optimized for use with multi-channel pipettes
- The refill system, depending on tip size, is packaged either dual-sided or in stack form
- Far less water compared to disposable boxes. All components 100% recyclable

Volume µl	Length mm	PU pc/Pack	Type epT.I.P.S. Reloads	epT.I.P.S. Reloads PCR clean
			Order No.	Order No.
0.1 – 10	34	960	10 69 73371	10 69 73754
0.1 – 20	40	960	10 69 73398	10 69 73770
0.5 – 20	46	960	10 69 73410	10 69 73797
2 – 200	53	960	10 69 73436	10 69 73819
20 – 300	55	960	10 69 73452	10 69 73835
50 – 1000	71	960	10 69 73479	10 69 73851
50 – 1250	76	960	10 69 73495	10 69 73878
50 – 1250	103	960	10 69 73630	10 69 73649
500 – 2500	115	480	10 69 73517	10 69 73894

❗ Available only in original packaging units. No fractions possible.

Pack: 10 trays of 96 tips (= 960 tips), or 10 trays of 48 tips (= 480 tips) for 500 – 2500 µl epT.I.P.S. Reloads



④ Pipette tips epTIPS Set

- epTIPS Set: 1 working box plus 5 refill trays
- Guard lid on refill trays enable contamination-free transfer of trays to the working box (epT.I.P.S. Box)
- Refill trays (epT.I.P.S. Reloads) can be autoclaved in advance
- System optimized for use with multi-channel pipettes
- The refill system, depending on tip size, is packaged either dual-sided or in stack form
- Far less waste compared to disposable boxes. All components 100% recyclable

Volume µl	Length mm	Order No.	PU pc/Set
0.1 – 10	34	10 69 73215	480
0.1 – 20	40	10 69 73231	480
0.5 – 20	46	10 69 73258	480
2 – 200 (2 – 100)	53	10 69 73274	480
20 – 300	55	10 69 73290	480
50 – 1000	71	10 69 73312	480
50 – 1250	76	10 69 73339	480
500 – 2500	115	10 69 73355	240

❗ Available only in original packaging units. No fractions possible.

Set: 1 box incl. 5 trays of 96 tips (= 480 tips), or incl. 5 trays of 48 tips (= 240 tips)



⑤ Pipette tips ep Dualfilter T.I.P.S.®

Eppendorf ep Dualfilter T.I.P.S.® are the first filter tips with contamination protection by means of a two-phase-filter. The distinctive, blue-white filter layers detain 100% of all aerosols and bio-molecules.

- Double protection by means of two filter layers
- Double protection for pipette and specimen
- Double protection from aerosols and bio-molecules
- Free of PCR inhibitor additives
- Sterile, Eppendorf PCR clean and pyrogen-free
- IvD-conform
- Charge certification

Version	Volume	Length mm	Order No.	PU pc/Pack
PCR	0.1 – 10 µl	34	10 69 78500	960
PCR	0.1 – 10 µl	40	10 69 78519	960
PCR	0.5 – 20 µl	46	10 69 78527	960
PCR	2 – 20 µl	53	10 69 78535	960
PCR	2 – 100 µl	53	10 69 78543	960
PCR	2 – 200 µl	55	10 69 78551	960
PCR	20 – 300 µl	55	10 69 78560	960
PCR	50 – 1000 µl	76	10 69 78578	960
PCR	0.25 – 2.5 ml	115	10 69 78586	240
PCR	50 – 1250 µl	103	10 69 78594	480
PCR	0.1 – 5 ml	120	10 69 78616	120
PCR	0.2 – 5 ml	175	10 69 78624	120
LoRentention PCR	0.1 – 10 µl	34	10 69 78632	960
LoRentention PCR	0.5 – 20 µl	46	10 69 78640	960
LoRentention PCR	2 – 100 µl	53	10 69 78659	960
LoRentention PCR	2 – 200 µl	55	10 69 78667	960
LoRentention PCR	20 – 300 µl	55	10 69 78675	960
LoRentention PCR	50 – 1000 µl	76	10 69 78683	960
Forensic	0.1 – 10 µl	40	10 69 78810	960
Forensic	2 – 20 µl	53	10 69 78829	960
Forensic	2 – 200 µl	55	10 69 78837	960
Forensic	50 – 1000 µl	76	10 69 78845	960
PCR	0.1 – 20 µl	42	10 69 78853	3840
PCR	2 – 100 µl	53	10 69 78861	3840

❗ Available only in original packaging units. No fractions possible.



⑥ Pipette tips epT.I.P.S. Racks Biopur

- Eppendorf Biopur pipette tips: Highest biological purity, e.g. guaranteed sterile, for RNA-se-, DNA-, ATP- and pyrogen-free work
- Each batch continuously controlled by external laboratory
- Batch-related quality certificate available over the Internet (www.eppendorf.com)
- Packed in boxes of 96, 48 or 24 tips (1 pack = 5 boxes)

Article	Volume	Length mm	Order No.	PU pc/Pack
epT.I.P.S.® Biopur, 0.1-20 µl	0.1 – 20 µl	40	10 69 75226	480
epT.I.P.S.® Biopur, 2-200 µl	2 – 200 µl	53	10 69 75234	480
epT.I.P.S.® Biopur, 20-300 µl	20 – 300 µl	55	10 69 75242	480
epT.I.P.S.® Biopur, 50-1000 µl	50 – 1000 µl	71	10 69 75250	480
epT.I.P.S.® Biopur, 50-1250 µl, 76 mm	50 – 1250 µl	76	10 69 75269	480
epT.I.P.S.® Biopur, 50-1250 µl, 103 mm	50 – 1250 µl	103	10 69 75277	480
epT.I.P.S.® Biopur, 0.25-2.5 ml	0.25 – 2.5 ml	115	10 69 75285	240
epT.I.P.S.® Biopur, 0.1-5 ml	0.1 – 5 ml	120	10 69 75293	120
epT.I.P.S.® Biopur, 0.2-5 ml	0.2 – 5 ml	175	10 69 75307	120
epT.I.P.S.® Biopur, 0.5-10 ml	0.5 – 10 ml	165	10 69 75315	120

❗ Available only in original packaging units. No fractions possible.



⑥

⑦ Pipette tips epT.I.P.S. Singles Biopur

- Eppendorf Biopur pipette tips, individually sealed in blister pack
- Batch number and use-by date on each blister pack
- Each batch continuously controlled by external laboratory
- Batch-related quality certificate available over the Internet (www.eppendorf.com)

Volume µl	Length mm	Order No.	PU pc/Pack
0.1 – 20	40	14 69 10019	100
2 – 200	53	14 69 10035	100
50 – 1000	71	14 69 10051	100

❗ Available only in original packaging units. No fractions possible.



⑦