

Graduated Pipettes, Class AS/Class B, BRAND

① Graduated pipettes, BLAUBRAND® ETERNA, class AS, type 3 (zero on top)

DE-M marked, with batch number and enclosed batch certificate. AR-Glas®. DIN EN ISO 835. Type 3 (zero on top), total delivery. Calibrated to deliver (Ex). Waiting time 5 seconds. Markings and inscriptions in ETERNA amber stain. On request, also available with individual certificate, USP individual certificate or DAkkS calibration certificate.

Volume ml	Graduation ml	Error limit ± ml	Order No.	PU pc/Pack
0.5	0.01	0.006	14 49 28421	12
1	0.01	0.007	14 49 28427	12
1 W	0.01	0.007	14 49 27427	12
2	0.02	0.010	14 49 28435	12
2 W	0.02	0.010	14 49 27435	12
5 W	0.05	0.030	14 49 28448	12
5 W	0.1	0.030	14 49 28449	12
10 W	0.1	0.05	14 49 28458	12
20 W	0.1	0.10	14 49 28463	6
25 W	0.1	0.10	14 49 28466	6

① Available only in original packaging units. No fractions possible.

W = with large suction tube suitable for cotton plugs

② Graduated pipettes, BLAUBRAND®, Class AS, type 2 (nominal volume on top)

DE-M marked, AR-Glas®. DIN EN ISO 835. Type 2 (nominal volume on top), total delivery. Calibrated to deliver (Ex). Waiting time 5 seconds. Markings and inscriptions in high contrast blue enamel.

All BLAUBRAND® graduated pipettes are supplied with batch number and one batch identification per packing unit of the manufacturer. On request, also available with individual certificate, USP individual certificate or DAkkS calibration certificate.

Volume ml	Graduation ml	Error limit ± ml	Overall length mm	Order No.	PU pc/Pack
0.5	0.01	0.006	360	10 49 27816	12
1	0.01	0.007	360	10 49 27817	12
1 W	0.01	0.007	360	10 49 27827	12
1	0.1	0.007	360	10 49 27818	12
2 ²⁾	0.01	0.010	360	10 49 27829	12
2	0.02	0.010	360	10 49 27819	12
2 W	0.02	0.010	360	10 49 27828	12
2	0.1	0.010	360	10 49 27820	12
5 W	0.05	0.030	360	10 49 27821	12
5 W	0.1	0.03	360	10 49 27822	12
10 W	0.1	0.05	360	10 49 27823	12
20 W	0.1	0.10	360	10 49 27824	6
25 W	0.1	0.1	450	10 49 27825	6
50 ¹⁾ W	0.5	0.2	450	10 49 27826	6

① Available only in original packaging units. No fractions possible.

¹⁾ in addition to ISO

²⁾ in addition to ISO size, not DE-M marked

W = with large suction tube suitable for cotton plugs

③ Graduated pipettes, BLAUBRAND® ETERNA, class AS, type 2 (nominal volume on top)

DE-M marked, AR-Glas®. DIN EN ISO 835. Type 2 (Nominal volume at the top), total delivery. Calibrated to deliver (TD, Ex). Waiting time 5 seconds. Ring mark and inscription in ETERNA amber stain.

All BLAUBRAND® graduated pipettes are supplied with a batch number and an accompanying batch certificate in the original packaging. On request, they are also available with an individual certificate, USP individual certificate or DAkkS calibration certificate.

Volume ml	Graduation ml	Error limit ± ml	Order No.	PU pc/Pack
5 W	0.1	0.03	10 49 27822	12
25 W	0.1	0.1	10 49 27825	6


① Available only in original packaging units. No fractions possible.

W = with large suction tube suitable for cotton plugs

Sale while stocks last.



WWQED052



- ④ **Graduated pipettes, BLAUBRAND®, USP, class AS, type 2 (nominal volume on top)**
DE-M marked, AR-Glas®. DIN EN ISO 835. Type 2 (nominal volumen on top), total delivery. Calibrated to deliver (Ex). Waiting time 5 seconds. Markings and inscriptions in high contrast blue enamel. Incl. USP batch certificate. On request, also available with USP individual certificate. Suitable for cotton plugs.

Volume ml	Graduation ml	Error limit ± ml	Overall length mm	Order No.	PU pc/Pack
1 ¹⁾	0.01	0.007	360	10 49 27506	12
2	0.02	0.01	360	10 49 27509	12
5	0.05	0.02	360	10 49 27511	12
10	0.1	0.03	360	10 49 27513	12

1) Error limit USP ± 0.01 ml

- ⑤ **Graduated pipettes, SILBERBRAND ETERNA, class B, type 2 (nominal volume on top)**
AR-Glas®. DIN EN ISO 835. Type 2 (nominal volumen on top), total delivery. Calibrated to deliver (Ex). Markings and inscriptions in ETERNA amber stain.

Volume ml	Graduation ml	Error limit ± ml	Order No.	PU pc/Pack
0.5	0.01	0.006	10 49 27816	12
1	0.01	0.010	10 49 27002	12
1	0.1	0.010	10 49 27004	12
2 ¹⁾	0.01	0.015	10 49 27005	12
2	0.02	0.015	10 49 27006	12
2	0.1	0.010	10 49 27820	12
5 W	0.05	0.05	10 49 27009	12
5 W	0.1	0.05	10 49 27010	12
10 W	0.1	0.08	10 49 27011	12
25 W	0.1	0.15	10 49 27013	6

1) in addition to ISO size, not DE-M marked
W = suitable for cotton plugs

