

The specialists: EM dispenser PP

Version with precision-ground glass piston for unproblematic aqueous solutions and for organic media that do not attack polypropylene, available in 2 popular sizes.

Operation

- „Click-stop-latch mechanism“ for reproducible variable volume settings.
- Sterilizable to max. 121°C

Conformity

- With Certificate of Conformity pursuant to DIN 12600
- Manufactured with state-of-the-art production engineering pursuant to DIN ISO 9001

Items supplied:

EM-dispenser PP, complete with screw thread A 32, 3 thread adapters (A 28, A 45 and S 40), 1 suction hose with union nut, 1 ejection unit, 1 valve mounting wrench, single calibration certificate.

Volume ml	Setting steps ml	Accuracy %	Variation coefficient ¹⁾ %	Order No.
2 – 10	0.25	0.6	0.2	14 37 34020
10 – 60	1.00	0.6	0.2	14 37 36020

¹⁾ Relative to the maximum volume.

