

Bottle-Top Dispenser VITLAB® piccolo, VITLAB

For dispensing tiny quantities of liquids in all areas of biochemical and medical research. Even the smallest quantities can be dispensed directly from the bottle with the VITLAB® piccolo – a big help, particularly for serial dispensing operations. Special advantage: Disposable tips are unnecessary. This reduces costs.

The ergonomic design makes dispensing effortless and stress-free. The VITLAB® piccolo can be operated with only one hand. Use the thumb to depress the volume dispensing button, just as with a pipette, and a reset mechanism refills the volume automatically. The discharge tube can be rotated over 360° so that it is always optimally situated with respect to the bottle label.

The bottle-top dispensers VITLAB®piccolo 1 and 2 are primarily used for aqueous and heavily diluted media. Checks are done in accordance with DIN EN ISO 8655-6 with the device filled to capacity and with uniform and jolt-free dosing.

- **VITLAB®piccolo 1** with a fixed volume
- **VITLAB®piccolo 2** with two fixed volumes

Included in delivery:

- ① VITLAB® piccolo 1 or 2 (GL 28 connecting threads), mounting tool, quality certificate and instruction manual.

Type	Volume µl	A ≤ ± %	CV ≤ %	Order No.
piccolo 1	100	3.0	0.4	10 85 31601
piccolo 1	200	2.5	0.4	10 85 31602
piccolo 1	250	2.0	0.4	10 85 31603
piccolo 1	500	1.5	0.3	10 85 31604
piccolo 1	1000	1.0	0.2	10 85 31605
piccolo 2	100/250	2.0	0.4	10 85 31606
piccolo 2	500/1000	1.0	0.2	10 85 31607
piccolo 2	1000/2000	1.0	0.2	10 85 31608

