Large Shakers, Bühler

Designed for heavy loads and continuous operation (floor placement only). Choice of orbital or reciprocal motion. Wide range of accessories for various applications. Includes basic unit without trays or stands. Please order the trays or stands that you need.

Technical data:

Max. load: 75 kg 30 – 200/min Frequency range: 760 × 600 mm Working surface:

Control: "control" with memory for shaking programs and touch-

programmable¹⁾ 1050 × 835 × 250 mm 100 – 240 V, 50/60 Hz Timer: External dimensions (W×D×H): Supply voltage:

VKS 75 Multi-flask shaker


	Article		Order No.
1	Multi-flask shaker VKS	S 75 A control	52 61 61110
	Shaking motion:	orbital	
	Shaking amplitude:	26 mm	

Article Order No. 2 Multi-flask shaker VKS 75 B control 52 61 61120

Shaking motion: reciprocal (horizontal)

Shaking amplitude: 50 mm or 36, 60, 80 mm (option)

	Tops and accessories for VKS 75 models			
	Article	Order No.		
	Oscillating plate VKS	52 61 52070		
	Base plate for fastening the VKS Combifix fittings or the VKS universal tray (1-level mode)			
	Combifix VKS Set A	52 61 51487		
	Frame rack with 9 clamping rails for upright vessels (1-level mode)			
3	Combifix VKS Set B	52 61 51488		
	Frame rack with 4 double clamping rails for tall or lying vessels (1 -level mode)			
	Combifix VKS Set C	52 61 51489		
	Frame rack with rails for fastening filter funnels (1-level mode)	32 01 31403		
	Universal tray VKS	52 61 51474		
	Tray with holes for spring clips or test tube racks (clips/racks on request) (1-level mode)			
	Top VKS 2-storey	10 61 62186		
	2-level frame rack for high loading capacities. Each level can be			
	loaded with either 1 tray size 760 × 600 mm or with 2 trays size			
	560×400 mm. Suitable trays on request. Distance between levels:			
_	400 mm. Height sufficient for 5000 ml Erlenmeyer flasks.			
4	Top VKS Gigant	52 61 52068		
	3-level frame rack for high loading capacities. Each level can be			
	loaded with either 1 tray size 760 × 600 mm or with 2 trays size			
	560 × 400 mm. Suitable trays on request. Distance between levels:			



340 mm. Height sufficient for 3000 ml Erlenmeyer flasks.









