

Blotting Chambers, JoJo Life Science



Electroblotting is a technique for immobilising proteins or for nucleic acid separation on a solid membrane support. The samples are then detected with specific antibodies, ligands or nucleic acid probes that bind individual proteins or nucleic acid sequences. We offer three types of systems: the modular electroblotters, tank transfer systems and semi-dry transfer systems.

omniBLOT blotter

- Mini-System transfers up to 4 Mini gels, Maxi-System up to 4 Maxi gels, 8 Mini Wide and 16 Mini gels
- Complete runtime flexibility from overnight to just a few hours
- Cool Pack for low-cost cooling for high-intensity transfer
- Compatible with OmniPage systems

omniBLOT blotter, high intensity


- Ideal for wet electroblotting of proteins - Western Blotting
- Platinum plate electrodes for increased transfer speed
- Fold-up cassettes for added convenience
- Suitable for gel strengths of 0.25 to 3 mm
- Compatible with OmniPage systems

electroBLOT blotter

- Ideal for wet electroblotting of proteins - Western Blotting
- Up to five gel-blot cassettes can be used in parallel
- Fold-up cassettes for added convenience
- Suitable for gel strengths of 0.25 to 3 mm

Semi-Dry blotter

- Fast transfer times
- Western-, Southern- and Northernblots
- Economical transfers thanks to very small buffer volumes
- Screw-on lid
- Gels of a thickness between 0.25 and 10 mm can be blotted
- Even heat distribution
- Electrodes with a long service life



Article	Order No.
① OmniBlot Mini, 10 × 10 cm blotting chambers incl. 4 cassettes	10 56 61519
② OmniBlot Mini Wide, 20 × 10 cm blotting chambers incl. 3 cassettes	10 56 61520
OmniBlot Maxi, 20 × 20 cm blotting chambers incl. 4 cassettes	10 56 61521
③ omniPAGE Blot Mini, 10 × 10 cm High intensity system including tank and lid, 2 cassettes, 8 fibre pads, cool pack	10 56 61523
WAVE Sub Blot, 20 × 20 cm High intensity system including tank and lid, 1 cassette, 6 fibre pads, cooling coil	10 56 61524
④ Electro Blot Mini, 10 × 10 cm Five Position incl. 5 cassettes	10 56 61522
Electro Blot Maxi, 20 × 20 cm Five Position incl. 5 cassettes	10 56 61091
⑤ Semi-Dry Mini, 10 × 10 cm System	10 56 61060
Semi-Dry Maxi, 20 × 20 cm System	10 56 60830

Please enter the webcode WW5RD001 to see the other variants in our online catalogue.

