

Micro Test Safe-Lock Eppendorf Tube®, Eppendorf



eppendorf

The patented hinged lid on these tubes eliminates any danger of leaks. Safe-Lock micro test tubes guarantee absolute safety during work with expensive or toxic samples, radioactive substances or DNA. The 0.5 ml, 1.5 ml and 2.0 ml test tubes have a graduation scale and a writing surface. The 0.5 ml Safe-Lock micro test tube is ideal for all PCR experiments.

- A small hook on the patented hinged lid clips around the rim of the test tube
- Tubes can be autoclaved when open (121°C, 20 min)
- Can be opened easily using one hand
- Frosted writing surface
- The lid hook prevents the tube from opening accidentally (e.g. during heating)
- Excellent mechanical stability during centrifugation
- Graduation scale
- Autoclavable micropestle enables resuspension of pellets in 1.5 – 2.0 ml test tubes
- Safe-Lock micro test tubes are also available in Eppendorf Biopur® and PCR clean quality



Colour	Version Eppendorf Quality™ Volume 500 pcs 0.5 ml Order No.	Version Eppendorf Quality™ Volume 1000 pcs 1.5 ml Order No.	Version Eppendorf Quality™ Volume 1000 pcs 2.0 ml Order No.
colorless	29 69 21023	29 69 20086	29 69 20094
yellow	29 69 21112	29 69 20159	29 69 20205
red	29 69 21120	29 69 20167	29 69 20213
blue	29 69 21139	29 69 20175	29 69 20221
green	29 69 21147	29 69 20183	29 69 20230
ambra (light protection)	29 69 21155	29 69 20191	29 69 20248
Color assortment	29 69 21708	29 69 21694	29 69 21686

❗ Available only in original packaging units. No fractions possible.

Assortment of yellow, red, blue, green and colorless.

